HiFil SupremeX

DE EN DA SV NO FI
PL LT ET NL CS

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Lichthärtbares NanoHybrid-Composite

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Light Curing NanoHybrid Composite

Product description

ORBIS HiFil SupremeX is a visible light curing radioactive nanohybrid composite for anterior and posterior restorations. It is based on Bis-GMA-resin and inorganic filler particles of 0.05 - 3 µm. The total filler load is 80%, the total filler volume is 69 %. Radiopacity of ORBIS HiFil SupremeX is 200 % aluminum (radiopacity of enamel: ~200 % aluminum, dentin: ~100 % aluminum). The composite meets the requirements of: DIN EN ISO 4049, type 1, class 2, group 1.

Indications/Intended use

Restorations of all cavity classes

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Application

Preparation of Tooth and Color Matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water or ORBIS Prophy Paste Vario Polish without Fluoride) prior to preparation and color matching.

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hellere Farben (z.B. A1; A2; C2) 20 Sekunden

dunkle Farben (z.B. A3; B3) 30 Sekunden

Hold the light emission window as close as possible to the filling material.

Finishing

Depending on the shape of the area contour with finishing diamonds, flexible grinding discs (e.g. ORBIPOL discs) or carbide burs, Polish to high gloss with polishing discs or polishing brushes (e.g. ORBIPOL micro-P, ORBIPOL Universal or ORBIPOL Spiral Polisher).

Storage

Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date. Close syringe immediately after use to avoid exposure to light.

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Storage

Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date. Close syringe immediately after use to avoid exposure to light.

Additional Notes/Warnings

- PLT are for single use only due to hygiene reasons.
- The ambient light of the dental lamp may start polymerization of the composite.
- Do not use any resin to adjust viscosity of the composite restorative.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.
- Keep away from children!

Composition

Dimethacrylates, dental glass, silicon dioxide, photo initiators

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUADMED - <https://ec.europa.eu/tools/eudamed>).

Warranty

First Scientific Dental GmbH guarantees that this product will be free from defects in material and manufacture. First Scientific Dental GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and First Scientific Dental GmbH's sole obligation shall be repair or replacement of the First Scientific Dental GmbH product.

Limitation

First Scientific Dental GmbH warrants this product will be free from defects in material and manufacture. First Scientific Dental GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and First Scientific Dental GmbH's sole obligation shall be repair or replacement of the First Scientific Dental GmbH product.

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