# Venus® Pearl

- GE** Gebrauchsanweisung
- EN** Instructions for use
- FR** Mode d'emploi
- ES** Instrucciones de uso
- IT** Istruzioni per l'uso
- PT** Instruções de uso

3609083071

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## Venus® Pearl

**Medical Device – for use by dental health care professionals only.**

### Intended purpose

Venus Pearl is a light curing, radio-opaque nano composite. It is indicated for adhesively bonded, tooth coloured restorations of anterior and posterior teeth including occlusal surfaces. Due to the excellent material properties, a beautiful aesthetic result can be easily reached. The radiopacity is equivalent to 2.0–2.5 mm of aluminum (dentin = 1 mm, enamel = 2 mm).

This dental material is suitable for patients requiring dental treatment for the following indications with consideration of the patient's medical history. Contact with the composite material during the dental restoration, for dental gaps and supports are due to the special situation of the treatment room in favour of the unsecured or infant.

Caution: Contact with the composite material during the dental restoration, for dental gaps and supports are due to the special situation of the treatment room in favour of the unsecured or infant.

Caution: Contact with the composite material during the dental restoration, for dental gaps and supports are due to the special situation of the treatment room in favour of the unsecured or infant.

Caution: Contact with the composite material during the dental restoration, for dental gaps and supports are due to the special situation of the treatment room in favour of the unsecured or infant.

**Composition**
Barium-Aluminium-Boro-Fluor-Silikate Glass, TCD-Urethancrylate, Silica, UDMA, TEGDMA, Polymer, Titanium Dioxid, Fluoreszenz Pigmente, Metalladditive, organische Pigmente, Ammoniumbenzoylperoxid, BHT, Kamphorquinon, Venus Pearl contains approximately 59% filler by volume (with 58% inorganic filler by volume), its filler particle size is 5 nm–5 µm.

### Clinical Indications

• Direct restoration of Class I–V cavities (according to G.V. Black)
• Direct composite veneers
• Shape correction of teeth (i.e. distal chatter, metal dependence, defects of the neck, etc.)
• Splinting of complete dentures, etc.
• Resin crowns
• Primary teeth
• Core build-up
• Repair of porcelain, composite (in combination with an adequate repair system)

### Contraindications

Application to dentine close to the pulp (less than 1 mm). The use of this product is contraindicated in case of known or suspected allergic components of this product. Do not use this product in case of known or suspected allergies against (meth)acrylate compounds.

### Potential side effects

This product or one of its components may in particular cases cause hypersensitive reactions.

### Warnings and safety instructions

Avoid by irritating to eyes and skin. Avoid contact with skin. May cause sensitization by skin contact. Wear suitable protective gloves. In case of contact with skin, immediately wash with plenty of water and soap. Avoid contact with eyes. In case of contact with eyes, rinse with water for several minutes. If symptoms persist, consult an eye specialist.

### Application hints

- Always inspect for damage before use. Damaged products must not be used. Composites should only be used at room temperature.
- In the case of deep cavity preparations in close proximity to the pulp or a direct pulp exposure, cover the dentine or the exposure locally with a calcium hydroxide liner covered with a glass ionomer cement before applying the adhesive.

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Do not use eugenol containing materials in or around the cavity preparation being restored, because these may interfere with the proper curing of the composite.							
Use adequate protection for spritzes to prevent from any contamination. Do not re-use after contamination.							
The PLTs are intended for single use only and reprocessing or re-use are not possible without damaging and contaminating the restorative material by those phenolic disinfection solutions. To avoid cross contamination PLTs must be disposed of immediately after initial use.							
When using Venus Pearl PLTs, use a gentle, even motion when exerting pressure on the applicator gun. Use of excessive force or sudden movement could result in unanticipated extrusion of the material or ejection of the PLT from the applicator gun, creating a potential hazard. We recommend using the Kulzer PLT Applicator STEEL. When using other PLT applicators gun see full instructions and precautions of these devices.							
<b>Available Shades</b>	1) Pure Shades for monochromatic layering						
	Corresponding VITA® shades						
<b>LIGHT</b>	A1, B1, B2, C1, D2						
<b>MEDIUM</b>	A2, A3, B3, C2, C3, C4, D3, D4						
<b>DARK</b>	A3.5, A4, B4						
<b>BLEACH</b>	BL, BL						
<b>ONE</b>	all 16 VITA® shades						
2) Art Shades for monochromatic and polychromatic layering							
<b>Opaque</b>	OB (Opaque Bleach)	OLC (Opaque Light Chromatic)	OMC (Opaque Medium Chromatic)	ODC (Opaque Dark Chromatic)	ODXC (Opaque Extra Dark Chromatic)		
	BL (Bleach Extra Light)	BL (Bleach Light)	A1	A2	A3	A3.5	A4
<b>Universal</b>		B1	C1	C2	C3	B2	B3
			D3			HKAZ.5	HKAS
<b>Incisal</b>	translucent	CL (Clear), AM (Amber)					
<b>opalescent</b>	CO (Clear Opal), YO (Yellow Opal)						
<b>Special</b>	core build-up	CORE (green)					
	gingiva	GUM (pink)					

Do not use eugenol containing materials in or around the cavity preparation being restored, because these may interfere with the proper curing of the composite.							
Use adequate protection for spritzes to prevent from any contamination. Do not re-use after contamination.							
The PLTs are intended for single use only and reprocessing or re-use are not possible without damaging and contaminating the restorative material by those phenolic disinfection solutions. To avoid cross contamination PLTs must be disposed of immediately after initial use.							
When using Venus Pearl PLTs, use a gentle, even motion when exerting pressure on the applicator gun. Use of excessive force or sudden movement could result in unanticipated extrusion of the material or ejection of the PLT from the applicator gun, creating a potential hazard. We recommend using the Kulzer PLT Applicator STEEL. When using other PLT applicators gun see full instructions and precautions of these devices.							
<b>Available Shades</b>	1) Pure Shades for monochromatic layering						
	Corresponding VITA® shades						
<b>LIGHT</b>	A1, B1, B2, C1, D2						
<b>MEDIUM</b>	A2, A3, B3, C2, C3, C4, D3, D4						
<b>DARK</b>	A3.5, A4, B4						
<b>BLEACH</b>	BL, BL						
<b>ONE</b>	all 16 VITA® shades						
2) Art Shades for monochromatic and polychromatic layering							
<b>Opaque</b>	OB (Opaque Bleach)	OLC (Opaque Light Chromatic)	OMC (Opaque Medium Chromatic)	ODC (Opaque Dark Chromatic)	ODXC (Opaque Extra Dark Chromatic)		
	BL (Bleach Extra Light)	BL (Bleach Light)	A1	A2	A3	A3.5	A4
<b>Universal</b>		B1	C1	C2	C3	B2	B3
			D3			HKAZ.5	HKAS
<b>Incisal</b>	translucent	CL (Clear), AM (Amber)					
<b>opalescent</b>	CO (Clear Opal), YO (Yellow Opal)						
<b>Special</b>	core build-up	CORE (green)					
	gingiva	GUM (pink)					

Do not use eugenol containing materials in or around the cavity preparation being restored, because these may interfere with the proper curing of the composite.					
Use adequate protection for spritzes to prevent from any contamination. Do not re-use after contamination.					
The PLTs are intended for single use only and reprocessing or re-use are not possible without damaging and contaminating the restorative material by those phenolic disinfection solutions. To avoid cross contamination PLTs must be disposed of immediately after initial use.					
When using Venus Pearl PLTs, use a gentle, even motion when exerting pressure on the applicator gun. Use of excessive force or sudden movement could result in unanticipated extrusion of the material or ejection of the PLT from the applicator gun, creating a potential hazard. We recommend using the Kulzer PLT Applicator STEEL. When using other PLT applicators gun see full instructions and precautions of these devices.					
<b>Available Shades</b>	1) Pure Shades for monochromatic layering				
	Corresponding VITA® shades				