ORBIS GI-Cem Luting

INSTRUCTIONS FOR USE
MD EU Medical device
In accordance with DIN EN ISO 9917-1

Application capsule

Product description:
ORBIS GI-Cem Luting is a radiopaque glass ionomer cement for luting in application capsules for direct application. ORBIS GI-Cem Luting is characterised by good biological properties, a high level of compressive strength, fluoride release and a low film-thickness.

Each capsule contains a minimum of 0.19 ml of useable cement.

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Application:
Carefully open the outlet by lifting the application tip until you feel a stop (see fig. 2). Place the opened application capsule into the holding bracket of an application instrument (see respective instructions for use). The cannula can be aligned as required by turning the whole application capsule in the holder. By operating the lever, dispense the material directly onto the surfaces which have been prepared for luting and distribute it using, for example, the application tip of the application capsule. Always apply a slight excess of material. Insert the restoration immediately applying light pressure and hold it in position continuing to apply light pressure while curing takes place. Avoid unnecessary waiting times between activation, mixing and application, since the incipient setting may make dispensing and the positioning of the restoration more difficult.

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Application:
At room temperature, ORBIS GI-Cem Luting allows a window of approx. 3 minutes working time, before setting begins after 3 - 5 minutes at 37 °C. Higher temperatures shorten and lower temperatures extend these working times. Prevent the ingress of water and saliva during application and setting phase. Excess material can be removed particularly easily during the rubber-elastic phase (before complete curing) using a scaler or a curette, for example.

Warnings, precautionary measures:
- Avoid direct application to dentine near the pulp or to exposed pulp.
- In case of eye contact, rinse well with plenty of water and consult an ophthalmologist.
- ORBIS GI-Cem Luting is not intended for non-reinforced all-ceramic restorations such as silicate ceramics, nor for indirect composite restorations, orthodontic brackets, ceramic inlays, ceramic onlays or veneers.

- Instruments used need to be rinsed immediately with cold water before residues of material cure. Following the removal of orthodontic bands, residues of material can be removed quickly and gently using an ultrasonic scaler.
- Ensure that the application capsule is mixed immediately after activation.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):
Strontium aluminium borosilicate glass, water, polyacrylylappo, Weinsäure

Storage instructions and application method:
Store ORBIS GI-Cem Luting in a dry place at temperatures of 15 °C to 23 °C. Only remove it from the foil immediately before use. Do not use after the expiry date.

Disposal:
Dispose of the product in accordance with local regulations.

Reporting obligation:
Serious events such as death, temporary or permanent serious deterioration of a patient's user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of ORBIS GI-Cem Luting must be reported to the manufacturer and the responsible authority.

Note:
The Summary of Safety and Clinical Performance of ORBIS GI-Cem Luting can be found in the European database on medical devices (EU-DAMED – <https://ec.europa.eu/tools/eudamed>).

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